

LSOP Title	Fluorescamine Method to Calculate Total Amino Acids
LSOP No.	LSOP12
Version	1.1
Location	UQ Node/Centre-wide
Policy/Procedure Link	
Risk Assessments	
Approved by	<i>Franziska Fichtner</i>
Date Approved	30/08/2021
Date Effective	03/06/2021
Next Review Date	03/06/2022
Contact for Assistance	plantsuccess@uq.edu.au

1.0 Scope

This procedure covers the protocol for calculating total amino acids using the Fluorescamine Method.

This CSOP does not cover the extraction from plant material.

2.0 Definitions

ddH₂O – double distilled water

RT – Room Temperature

Na-Borate – Sodium Borate

3.0 Material and Equipment

1. Synergy Microplate Reader
2. Glutamate
3. Fluorescamine
4. ddH₂O
5. Na-Borate



6. Pipette & Pipette Tips

7. Eppendorf Tubes



4.0 Prescribed Actions

1. Create standards with glutamate at concentrations of

- a. 0 uM
- b. 100 uM
- c. 200 uM
- d. 400 uM
- e. 600 uM
- f. 800 uM

2. Duplicate 4 uL aliquots of the filtered aqueous extract from the methanol/chloroform extraction.



3. Add to standards and duplicates:

- a. 100 uL ddH₂O
- b. 15 uL 0.1 M Na-Borate (pH 8)
- c. 90 uL of 0.1% (w/v) fluorescamine (in acetonitrile)

4. Incubate at RT for 5 min

5. Fluorescence at 485 nm in a Synergy Microplate Reader with excitation at 405 nm.

6. Use a standard curve to calculate your total amino acids in nmol/gFW