

LSOP Title	Tre6P Sugar Extraction
LSOP No.	LSOP34
Version	1.1
Location	UQ Node/Centre-wide
Policy/Procedure Link	UQ- Equipment OGTR
Risk Assessments	
Approved by	<i>Franziska Fichtner</i>
Date Approved	30/08/2021
Date Effective	30/06/2021
Next Review Date	30/06/2022
Contact for Assistance	plantsuccess@uq.edu.au

1.0 Scope

To outline the procedures for Chloroform/Methanol extractions of all water-soluble metabolites (including trehalose 6-phosphate)

2.0 Definitions



Eppendorf tube –

Tre6P – Trehalose-6-Phosphate

3.0 Materials and Equipment

1. Speed vac
2. Vortex
3. Eppendorf tubes
4. Micro pestle
5. Ice
6. Liquid nitrogen
7. Ground plant material



8. Weigh balance
9. Chloroform/Methanol
10. Pipette and Pipette Tips
11. Ice cold water



4.0 Prescribed Actions

Chloroform/Methanol extraction

1. weigh ~18 mg (15 to 20 mg range) ground plant material in 2 ml safe-lock Eppendorf tubes (we can also use single pea buds, then the elution volume at the end changes)
2. Add 350 μ l Chloroform/Methanol (3/7; icecold, store at -20°C) into the tube which still is placed in liquid nitrogen and mix well.
3. After having done 10 samples vortex all and put into -20°C . Continue with the rest of the samples.
4. After 2 h in -20°C , place the samples on ice again and add 350 μ l of ice-cold water, vortex and centrifuge for 10 min at 4°C .
5. Transfer the upper aqueous phase into a new tube (1 ml screw cap)
6. Re-extract the lower Chloroform phase with 300 μ l of ice-cold water.
7. Combine the water phases and evaporate to dryness in the speed vac.
8. Dissolve samples in 400 μ l water (this depends on how much material we have and can change, for single buds use less)



For Tre6P measurements on LC-MS/MS

- filter through a MultiScreen PCR-96 Filter Plate membranes (Merck Millipore, www.merckmillipore.com) to remove high molecular weight contaminants.