

LSOP Title	Using the plate reader
LSOP No.	LSOP49
Version	1.1
Location	UQ Node/Centre-wide
Policy/Procedure Link	<a href="#">UQ- Equipment</a>
Risk Assessments	
Approved by	<i>Lindsay Shaw</i>
Date Approved	13/10/2021
Date Effective	05/08/2021
Next Review Date	05/08/2022
Contact for Assistance	<a href="mailto:plantsuccess@uq.edu.au">plantsuccess@uq.edu.au</a>

## 1.0 Scope

To outline the use of the plate reader.

## 2.0 Prescribed Actions

- Turn plate reader and PC on (Fig 1)




Figure 1: Power button for the plate reader and PC

**LABORATORY STANDARD OPERATING PROCEDURE (LSOP)**

ARC COE for Plant Success in Nature and Agriculture: *Using the Plate Reader*

*NB: Plate reader will make some sounds while it calibrates*

- Plate reader is ready when the black tray is PERMANENTLY located outside (see Fig 2). It will move in and out while calibrating.



Figure 2: Plate reader with the black tray out

- Start Gen5 software (see button circled in Fig. 3) on the computer and select your protocol of choice




Figure 3: Gen5 software

- Put plate on the tray and follow instructions in software to start reading.
- To shut down first press the black button on the reader (red circle Fig 4.) and wait for black tray to go inside the reader and then press the off button (see blue circle Fig 4.)

**LABORATORY STANDARD OPERATING PROCEDURE (LSOP)**

ARC COE for Plant Success in Nature and Agriculture: *Using the Plate Reader*



Figure 4: Tray (red circle) and OFF (blue circle) button on the plate reader

- Shut down PC