

LSOP Title	Arabidopsis growth on agar plates
LSOP No.	LSOP03
Version	1.1
Location	UQ Node/Centre-wide
Policy/Procedure Link	<a href="#">UQ- Equipment</a> <a href="#">OGTR</a>
Risk Assessments	
Approved by	Lindsay Shaw
Date Approved	12/10/2021
Date Effective	02/06/2021
Next Review Date	02/06/2026
Contact for Assistance	<a href="mailto:plantsuccess@uq.edu.au">plantsuccess@uq.edu.au</a>

## 1. Scope

*To outline the procedures for growing Arabidopsis on agar for CK extraction*

*This LSOP does not cover the growth of Arabidopsis in soil or any other medium. It also does not include how to create agar or agar plates.*

## 2. Definitions



EtOH – ethanol

Murashige and Skoog (MS) – a type of selective medium

## 3. Materials and Equipment

1. Square petri dishes
2. MS agar
3. Sucrose
4. Ethanol
5. Pipette
6. Microcentrifuge tube
7. Distilled water
8. Alfoil



## 9. Laminar Flow

## 4. Prescribed Actions



1. Prepare square petri dishes with 50 mL half-strength Murashige and Skoog (MS) medium containing 1% sucrose and 1% agar (pH 5.7).

*NB: autoclave 20% sucrose solution and MS agar separately, then add sucrose to the MS medium prior to pouring the plates.*

*NB: you can create any media you require for this step*

2. Sterilise *Arabidopsis* seeds (LSOP05)
3. Dispense seeds to be used into a new sterile microcentrifuge tube and add ~200  $\mu$ L of sterile distilled water or 2% agar (see appendix or LSOP05 appendix)
4. Using a pipette, suck up water/agar until you get a single seed in the pipette tip, then dispense the liquid + seed onto the plate. Place seeds in a line across the center of the plate. Repeat until you have dispensed ~15-20 seeds per plate across all plates.
5. Wrap plates in alfoil and stratify by keeping at 4°C for two days.  
*N.B. This should be done with plates flat.*
6. Places plates upright (with the line of seeds horizontal) in a growth cabinet and grow for 4 weeks at desired growth conditions (i.e. 150, 25°C, 12h days).

## 5. Appendix

## 0.2% Agar

1. In a small Schott bottle, add 0.2 g agar to 100 mL milliQ water.
2. Autoclave (LSOP83).
3. Add ~ 1.5 mL (depending on number of seeds) sterile 0.2% agar solution to *Arabidopsis* seeds in Eppendorf to vernalise / plate