

LSOP Title	Arabidopsis Seed Surface Sterilisation
LSOP No.	LSOP05
Version	1.1
Location	UQ Node/Centre-wide
Policy/Procedure Link	<a href="#">UQ- Equipment</a> <a href="#">OGTR</a>
Risk Assessments	
Approved by	Franziska Fichtner
Date Approved	30/08/2021
Date Effective	28/05/2021
Next Review Date	28/05/2026
Contact for Assistance	<a href="mailto:plantsuccess@uq.edu.au">plantsuccess@uq.edu.au</a>

## 1.0 Scope

*This procedure covers the sterilisation protocol for Arabidopsis seeds surfaces to reduce contamination and growth of bacteria.*

*This LSOP does not cover Arabidopsis growth methods or protocols, seed removal from mature Arabidopsis plants or anything unrelated to seed surface sterilisation.*

## 2.0 Definitions

*Arabidopsis – Arabidopsis thaliana*



EtOH – Ethanol

milliQ water – high purity water from the milliQ system

Tween20 – a specific type of polysorbate 20 surfactant (i.e. detergent)

### 3.0 Materials and Equipment

1. Eppendorf tube

2. Falcon tube

3. Pencil

4. Bleach

5. Distilled water

6. Tween20

7. Ethanol

8. Pipette

9. Centrifuge

10. Laminar Flow



### 4.0 Prescribed Actions

1. Decant *Arabidopsis* seeds into Eppendorf tubes, labeled with pencil (doesn't wash off in ethanol)

2. Make sterilizing solution (must be made fresh each time) in a falcon tube as follows:

a. 25 mL bleach (12.5% strength; kept under sink)

b. 25 mL distilled water

c. Drop of Tween20 (dispense drop using a 200uL pipette)

d. Invert to mix.



3. In the laminar flow (or lab), add 1 mL 70% EtOH to Eppendorf tube with *Arabidopsis* seeds. Wash down ethanol is fine to use. Invert tube a few times to rinse the seeds, then quickly discard EtOH (seeds in ethanol too long will start to die).



4. Add 1 mL of sterilizing solution to seeds in Eppendorf tube. Invert to mix at least four times over a 10 min timed period (period may be extended up to 15 minutes if there are greater than 500uL of seeds or if the seeds came from plants with a lot of dirt or fungus, but do not exceed 10-15 minutes as seeds will die).



5. In laminar flow, remove sterilizing solution, then rinse seeds with 1 mL sterile distilled water (autoclaved) at least four times, using fresh water each rinse, and at least until detergent bubbles are not visible.

*NB: Can use a brief centrifuge to help stop seeds from floating during rinses if required.*

6. Add 0.2% agar (sterile; autoclaved; see appendix) to *Arabidopsis* seeds in Eppendorf tubes, brief centrifuge to remove any air bubbles if needed.

## 5.0 Appendix

### 0.2% Agar

1. In a small Schott bottle, add 0.2 g agar to 100 mL milliQ water.
2. Autoclave.
3. Add ~ 1.5 mL (depending on number of seeds) sterile 0.2% agar solution to *Arabidopsis* seeds in Eppendorf to vernalise / plate