

LSOP Title	DNA extraction (quick and dirty)
LSOP No.	LSOP08
Version	1.1
Location	UQ Node/Centre-wide
Policy/Procedure Link	UQ- Equipment UQ -waste OGTR
Risk Assessments	
Approved by	Francois Barbier
Date Approved	14/09/2021
Date Effective	03/06/2021
Next Review Date	03/06/2026
Contact for Assistance	plantsuccess@uq.edu.au

1.0 Scope

This protocol outlines the procedures for completing DNA extraction on leaf tissue.

2.0 Definitions

RT – Room Temperature

PCR – Polymerase Chain Reaction



TRIS – Trisaminomethane



NaCl – Sodium Chloride



EDTA – Ethylenediaminetetraacetic acid



SDS - Sodium Dodecyl Sulfate

NB: can also be known as SLS or Sodium Laurel Sulfate

3.0 Materials and Equipment

1. Plant tissue
2. Extraction Buffer (see appendix)
3. Centrifuge
4. Isopropanol





5. Ethanol
6. Water
7. Eppendorf tubes
8. Pipette (& Pipette Tips)

4.0 Prescribed Actions

1. Harvest 1 leaf per plant

NB: can be stored in -20°C for months



2. Add 200 µL of Extraction Buffer (see appendix) to your sample and grind (e.g. using a tip or tooth pick)

3. Spin at 10 min at full-speed and RT



4. Take 150 µL of the supernatant and mix with 150 µL isopropanol

5. Incubate for 10 min at RT

NB: here you can also take a break and put samples in the fridge

6. Spin 10 min at full-speed and RT



7. Discard supernatant and wash with 70% ethanol (v/v)
8. Dry pellet and resuspend in 80 µL water
9. Take 2 µL for a PCR reaction (you may have to dilute if the sample is too concentrated)

5.0 Appendix

Extraction Buffer (EB):



200 mM TRIS (pH 7.5)



250 mM NaCl



25 mM EDTA (pH 8.0)



0.5% SDS