

LSOP Title	Soil screening of primary transformants (only for BASTA resistance)
LSOP No.	LSOP29
Version	1.1
Location	UQ Node/Centre-wide
Policy/Procedure Link	<a href="#">UQ- Equipment</a> <a href="#">UQ -waste</a> <a href="#">OGTR</a>
Risk Assessments	
Approved by	Milos Tanurdzic
Date Approved	11/10/2021
Date Effective	18/06/2021
Next Review Date	18/06/2026
Contact for Assistance	<a href="mailto:plantsuccess@uq.edu.au">plantsuccess@uq.edu.au</a>

## 1.0 Scope

*This procedure covers the soil screening of primary transformants (only for BASTA resistance) based on the floral dipping of Arabidopsis based on the Zhang et al., 2006 short version method*

*This LSOP does not cover floral dipping of other species.*

## 2.0 Definitions

Cotyledons – embryonic leaf

## 3.0 Materials and Equipment

1. Spray bottle
2. BASTA or glufosinate ammonium
3. Clear plastic film

## 4.0 Prescribed Actions

1. Spread non-sterile seed onto moistened soil under normal condition



## LABORATORY STANDARD OPERATING PROCEDURE (LSOP)

ARC COE for Plant Success in Nature and Agriculture: *Soil screening of primary transformants (only for BASTA resistance)*

2. Cover the soil and tray with a clear plastic film until the germinated seedlings develop four to six leaves
3. Initiate the herbicide spraying when the cotyledons are visible, normally around 8-10 d after sowing.
4. Spray the seedlings two times per week with a diluted solution of 250 mg L<sup>-1</sup> herbicide (BASTA or glufosinate ammonium).



*NB: True transformants will develop, while non-transformants will become chlorotic and eventually die after 3-4 weeks of herbicide treatment.*

## 5.0 Appendix

Read the full protocol before starting with this short version (the side notes are useful)

Zhang, X., et al. (2006) Agrobacterium-mediated transformation of *Arabidopsis thaliana* using the floral dip method.

<https://www.nature.com/articles/nprot.2006.97>